



UNICORN LYTE

UNIQUE ELECTROLYTE ANALYZER

UNICORN LYTE 3 (Na⁺ K⁺ Cl⁻)



UNICORN LYTE 5 (Na⁺ K⁺ Cl⁻ Ca²⁺ pH)



UNIQUE TECHNOLOGY

- Advance Multi Biosensor Technology (MBT)
- No maintenance for electrodes

EASE OF USE

- Closed Reagent Pack easy to change
- Low sample volume

HIGHER STABILITY

- Durable Ceramic Rotary Valve
- Concise Liquid Flow Path
- Various auto Calibration

UNICORN LYTE

UNIQUE ELECTROLYTE ANALYZER

Specification

Method | Multi-Biosensor Technology



Sample | 60-120µl

Speed | ISE≥80 tests/hour

Measure scope & Precision

Parameters	Detection Range (mmol/L)	Precision (CV)	Resolution (mmol/L)
K ⁺	0.50 ~ 20.00	≤1.0%	0.01
Na ⁺	15.0 ~ 200.0	≤1.0%	0.1
Cl ⁻	15.0 ~ 200.0	≤1.0%	0.1
Ca ⁺	0.10 ~ 200.0	≤1.0%	0.01
pH	4.00 ~ 9.00	≤1.0%	0.01

Printer	Inner thermal printer
Expand Function	RS232 & COM Connector, Computer and Barcode reader available
Ambient Conditions	0-40 °C, Relative humidity ≤ 85%
Power	AC 80V ~ 260V, 50-60Hz, ≤ 80W
Power Requirement	Earth Voltage <5V, Online UPS (230V) with isolated, transformers - 0.5 KVA



TOLL FREE
1800-270-1027



EC REP

UNIONPROJET S.R.L.S.
Via Aliprandi 41-20851
Lissone (MB) Italy

ISO 9001 & 13485
Certified Company



Exclusively Marketed in India by:

VECTOR BIOTEK PVT. LTD.

Factory & Registered Office:
425, NEW GIDC, KABILPORE, NAVSARI - 396 424. (GUJ.) INDIA
Toll free:1800-270-1027 | Tel : +91-93757 94103 | Fax: +91-2637-265105
E-mail : vector.support@beaconindia.com | Web: www.vectorbiotekindia.com
CIN No. U33112GJ2002PTC040978