

BPL Medical
Technologies



Experience the Trace Quality

THREE CHANNEL ECG RECORDER
CARDIART 6208 View



Happier Living Everyday

www.bplmedicaltechnologies.com

CE 0843



Enhanced ECG diagnosis made possible with superior trace quality



Three Channel ECG Recording

Simultaneous 3 channel ECG recording with 12 lead simultaneous acquisition and Arrhythmia Detection in auto mode.



Large LCD Display with Patient Data Entry feature

Large Colour LCD Display enabling the preview of ECG and preview of patient information entry.



Internal Memory

Internal memory which allows storage of up to 200 ECGs.



Ergonomic Design

The compact, light weight (only 2.5 kg) ensures device portability



Power Line Safe Design

Designed for safe use, both on power line as well as battery



Built-in Rechargeable Battery

High performance Li-ion rechargeable battery enables mains independent operations and provides battery use after charging



Digital Filtering

Removes interference from power line, muscle tremor, etc. with automatic base line correction



USB

Enables the storage and transfer of ECG recordings from an ECG machine to PC



Inbuilt PDF converter

Facilitates transfer of ECGs from the device to USB in PDF format



Optional

Data transfer to PC through Ethernet port, Pacemaker detection

CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
CIN: U33110KA2012PTC067282
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in



A BPL Medical Technologies Magazine



Follow Us on

